

## ASTAR - PhysioMG 827

### Magnetoterapia



#### ERGONOMY

- treatments in Dual Mode
- visual indicator of the magnetic field activity – illuminator
- a sensory indicator of the magnetic field activity – magnetic stripes
- setting the frequency range of the magnetic field with an accuracy of 0,01 Hz using the on-screen keyboard
- setting of the time value of the treatment with an accuracy up to 1 sec using the on-screen keyboard
- two independent treatment channels
- 7" color touchscreen
- operation mode: program/ manual
- built-in treatment programs
- disease entities selected by name or field
- illustrated encyclopedia with a description of treatment methodology
- a handy list of favorite programs
- the ability to edit program names and user sequences
- statistics of the performed treatments
- volume control of the signaling sound
- self-test – current check of the unit performance

#### MAGNETIC FIELD

- field shape: rectangular modified, triangular, sinusoidal, impulse, half rectangular modified, semi-triangular, half-sine, half-pulse
- continuous and impulse emission
- wide frequency range
- high induction of the magnetic field

#### MANUAL MODE

full control over the treatment parameters for the advanced users

#### TREATMENT PROGRAMS

SPORTS PHYSIO S.r.l. sede legale C.da Torrazze, snc – 95127 Catania

- simplified operation on the unit
- 50 built-in treatment programs for the CS75 applicator
- 50 built-in treatment programs for the CS60 applicator
- 28 built-in treatment programs for the CS35 applicator
- 46 built-in treatment programs for the CP applicator
- 45 built-in treatment programs for the CPEP applicator
- 44 built-in treatment programs for the dual mode
- 50 user-defined programs
- favorite programs

## TREATMENT SEQUENCES

- simplified operation on the unit
- 10 user-defined sequences
- creating your own names of the sequences

## TECHNICAL PARAMETERS

- maximum nominal induction in the geometric center of the applicator:
  - applicator CS75 3 mT
  - applicator CS60 4 mT
  - applicator CS35 10 mT
  - applicator CP 3 mT
  - applicator CPEP 25 mT
- maximum nominal induction at the surface of the applicator:
  - applicator CS75 4,5 mT
  - applicator CS60 6 mT
  - applicator CS35 12 mT
  - applicator CP 12,5 mT
  - applicator CPEP 50 mT
- the maximum value of induction change:
  - applicator CS75 9 mT
  - applicator CS60 12 mT
  - applicator CS35 24 mT
  - applicator CP 25 mT
  - applicator CPEP 100 mT
- parameters of the intermittent mode pulse 1 s/ pause 0,5 – 8 s
- frequency spectrum 0 – 50 Hz
- treatment clock maximum 60 minutes
- power supply, power consumption 230 V  $\pm$ 10 %, 50/60 Hz, 350 VA
- controller dimensions 34 x 28 x 11 cm
- controller weight max 7 kg

## STANDARD EQUIPMENT

- Mains cable 1pc
- Spare fuses WTA-T 3,15 A/250 V 2pcs
- Permanent magnet 1pc
- Protective glasses for the patient – goggles 1pc
- Tracer for a resistive touch screen 1pc
- Cloth for the LCD screen 1pc
- Magnetic stripe 1pc


**ACCESSORI DISPONIBILI**



Applicatori CPEP con fasce elastiche



Solenoida da 35 cm - CS35



Solenoida da 60 cm - CS60;



Solenoida da 75 cm – CS75;



Letto per Magnetoterapia



Applicatori CP con carrello



---

SPORTS PHYSIO S.r.l. sede legale C.da Torrazze, snc – 95127 Catania



Carrello VERSA



Carrello per solenoide CS35

