

ASTAR – PhysioGo 200 A

Ultrasuono



ERGONOMICS

- operation mode: program / manual
- preset treatment program base
- statistics – the number of performed treatments
- encyclopedia with treatment methodology description
- treatment programs chosen by name or medical field
- concise list of favorite programs
- buzzer volume setting
- possibility of edition the names of user-defined programs
- calibration parameters stored in the ultrasound head
- self-test – systematic control of the unit
- 7" color touchscreen
- real time clock

ULTRASOUND THERAPY

- operation frequency 1 MHz/3,5 MHz
- effective radiation area 1 cm² and 4 cm²
- continuous emission – thermal effect
- pulse emission – micro massage effect
- ultrasound heads provided with the contact control connected to treatment timer
- precise control over the ultrasound wave emission during performed treatment
- head sensitivity calibration
- automatic control of emitted power

MANUAL MODE

Complete control over parameters for advanced users

TREATMENT PROGRAMS

- 58 preset treatment programs
- 50 user-defined programs
- 58 favorite programs



TECHNICAL SPECIFICATION

- | | |
|--|---|
| <ul style="list-style-type: none"> ultrasound generator operation mode continuous, pulse max ultrasound wave intensity in continuous / pulse mode pulse frequency duty factor for pulse mode mains supply dimensions weight | <p>2 / 3 W/cm²
 16 Hz, 48 Hz, 100 Hz
 5 – 75 %, step 5 %
 230 V, 50/60 Hz, max 75 W, 90 VA
 34 x 28 x 11-16 cm
 6 kg</p> |
|--|---|

STANDARD ACCESSORIES

- ultrasound head 1/ 3,5 MHz, 4 cm² with holder
- mains cable
- gel 500 g (1 pc)
- spare fuses WTA-T1 A/250 V (2 pcs)
- touchscreen stylus pen (1 pcs)
- touchscreen cleaning cloth (1 pcs)
- masking covers with cutout (2 pcs)

ACCESSORI DISPONIBILI

<input type="checkbox"/> Ultrasound head 1/ 3,5 MHz, 1 cm ² with holder 	<input type="checkbox"/> Carrello VERSA 
--	---